


RaktDrishti (TM) Calibration Reference Card


OMNIKON 2026 BIOTECH - NON-INVASIVE ANEMIA SCREENING AID



Frontline Health Worker Instructions (ASHA / ANM / Anganwadi):

1. Position this reference card beside the eye conjunctiva, fingernails, or palm.
2. Keep all color patches and corner markers visible and unshadowed in camera frame.
3. Hold device ~15-20 cm away under indirect ambient lighting. Avoid flash glare.

IMPORTANT: Print at 100% scale (Do NOT fit-to-page). Validated screening aid.



CARD-ID: RD-CALIB-2026-V1.0 | MATRIX: D65-sRGB | SCALE: 100%